

# Anti-HEXA (1G3) Mouse antibody

Catalog Number: 168466

## Product name

Anti-HEXA (1G3) Mouse antibody

## Specificity

Human

## Antibody description

HEXA (1G3) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human HEXA expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications

WB, FC


## Dilutions

WB: 1/500 - 1/2000


FC: 1/200 - 1/400

ELISA: 1/10000


## Validations



Western blot detection of HEXA in Jurkat, MCF7 and C6 cell lysates using HEXA mouse mAb (1:2000 diluted). Predicted band size: 61 kDa. Observed band size: 50 kDa.



Western blot analysis using HEXA mouse mAb against L1210 (1), and HL7702 (2) cell lysate.



Flow cytometric analysis of HeLa cells using HEXA mouse mAb (green) and negative control (red).