

PHYHIP antibody



Catalog Number: 113866

Product name

PHYHIP antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PHYHIP Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 35 kDa, 43-45 kDa

Preparation

This antibody was obtained by immunization of PHYHIP recombinant protein (Accession Number: XM_006716416). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113866(PHYHIP antibody) at dilution of 1:1000