

Casein Kinase 1 epsilon antibody



Catalog Number: 108866

Product name

Casein Kinase 1 epsilon antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Casein Kinase 1 epsilon Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human heart tissue, human kidney tissue, human ovary tissue, human testis tissue.

Preparation

This antibody was obtained by immunization of Casein Kinase 1 epsilon recombinant protein (Accession Number: NM_152221). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:108866(CSNK1E antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:108866(CSNK1E antibody) at dilution of 1:50 (under 40x lens)