

# Human FEN1/RAD2 recombinant protein



Catalog Number: 519866

## General Information

### Gene Name Synonym

Flap Endonuclease 1, FEN-1, DNase IV, Flap Structure-Specific Endonuclease 1, Maturation Factor 1, MF1, hFEN-1, FEN1, RAD2

### Protein Construction

A DNA sequence encoding the Human FEN1/RAD2 Met1-Lys380 is expressed.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

## Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

## Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 50 mM Tris, 50mM NaCl, 1mM DDT, 10% Glycerol, pH 8.0.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

## SDS-PAGE



Human FEN1/RAD2 recombinant protein