

# TRPV3 Polyclonal Antibody



Catalog Number: 161166

## Product name

TRPV3 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to TRPV3

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 400-500 of human TRPV3 (NP\_659505.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

## Validations



Immunohistochemistry - TRPV3 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat heart using TRPV3 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - TRPV3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using TRPV3 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - TRPV3 Polyclonal Antibody

# TRPV3 Polyclonal Antibody



Catalog Number: 161166

---

Immunohistochemistry of paraffin-embedded

mouse brain using TRPV3 antibody at dilution of 1:100 (40x lens).