

# ZNF354A antibody



Catalog Number: 117166

## Product name

ZNF354A antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ZNF354A Rabbit Polyclonal antibody. Positive WB detected in DU 145 cells. Observed molecular weight by Western-blot: 69 kDa

## Preparation

This antibody was obtained by immunization of ZNF354A recombinant protein (Accession Number: XM\_047417644). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



DU 145 cells were subjected to SDS PAGE followed by western blot with Catalog No:117166(ZNF354A antibody) at dilution of 1:1000