

MCUR1 antibody

Catalog Number: 112566

Product name

MCUR1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

MCUR1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HEK-293 cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 40-50 kDa

Preparation

This antibody was obtained by immunization of MCUR1 recombinant protein (Accession Number: NM_001031713). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112566(CCDC90A Antibody) at dilution of



1:1000

IP Result of anti-CCDC90A (IP:Catalog No:112566, 3ug; Detection:Catalog No:112566 1:300) with HEK-293 cells lysate 920ug.