

Anti-Calcineurin A Rabbit antibody



Catalog Number: 167566

Product name

Rabbit IgG

Anti-Calcineurin A Rabbit antibody

Applications

Specificity

WB, FC, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/2000

Calcineurin A Rabbit monoclonal antibody

FC: 1/20

Preparation

IP: 1/20

Antigen: A synthetic peptide of human Calcineurin A

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of Calcineurin A in HeLa, A549, U2OS cell lysates using Calcineurin A Rabbit mAb(1:1000 diluted). Predicted band size: 59kDa. Observed band size: 59kDa.