

Human UBE2J2/NCUBE2 (GST tag) recombinant protein



Catalog Number: 518566

General Information

Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 J2, Non-Canonical Ubiquitin-Conjugating Enzyme 2, NCUBE-2, UBE2J2, NCUBE2

Protein Construction

A DNA sequence encoding the Human UBE2J2/NCUBE2 Glu124-Arg224 is expressed with a GST tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.2.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2J2/NCUBE2 (GST tag) recombinant protein