

Human ANGPTL4 / angiopoietin-like protein 4 (His Tag)



Catalog Number: 503966

General Information

Gene Name Synonym

Angiopoietin-like protein 4; Hepatic fibrinogen/angiopoietin-related protein

Protein Construction

A DNA sequence encoding the human ANGPTL4 (Q9BY76) (Pro166-Ser406) was expressed with an N-terminal polyhistidine tag.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human ANGPTL4 comprises 261 amino acids and has a predicted molecular mass of 29.5 kDa. The apparent molecular mass of the protein is approximately 35 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human ANGPTL4 Protein (His Tag) SDS-PAGE