

Anti-ABAT/GABA-T Rabbit antibody



Catalog Number: 166966

Product name

Anti-ABAT/GABA-T Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ABAT/GABA-T Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human GABA-T

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/100

ICC/IF: 1/20

FC: 1/100

IP: 1/20

Validations



Western blot detection of GABA-T in Rat Brain lysates using GABA-T Rabbit pAb(1:1000 diluted). Predicted band size:56kDa. Observed band size:56kDa.