

Anti-phospho-Rb (Ser608) Rabbit antibody



Catalog Number: 167966

Product name

Anti-phospho-Rb (Ser608) Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, IHC-P

Antibody description

phospho-Rb (Ser608) Rabbit monoclonal antibody

Dilutions

WB: 1/2000-1/10000

Preparation


Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser608 of human Rb

IHC: 1/20

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Validations



Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Rb (phospho-Ser608) in C6,Ramos using Rb (phospho-Ser608) Rabbit mAb(1:1000 diluted)

Clonality

Monoclonal