

Anti-PAX6 Rabbit antibody



Catalog Number: 168966

Product name

Rabbit IgG

Anti-PAX6 Rabbit antibody

Applications

Specificity

WB, IHC-P, FC

Human, Mouse, Rat

Dilutions

Antibody description

1:1000-1:5000

PAX6 Rabbit monoclonal antibody

1:20-1:100

Preparation


1:100-1:500

Antigen: A synthetic peptide of human PAX6

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of PAX6 in Rat Brain lysates using PAX6 Rabbit mAb(1:1000 diluted).Predicted band size:47kDa.Observed band size:47kDa.

Clonality

Monoclonal

Ig Type