

Anti-POTEB3 antibody



Catalog Number: 176266

Product name

Anti-POTEB3 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to POTEB3

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human a26b1. aa range:481-530

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western blot analysis of Mouse-kidney 3T3 lysis using POTE15 antibody. Antibody was diluted at 1:500