

TPK1 antibody

Catalog Number: 116266

Product name

TPK1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TPK1 Rabbit Polyclonal antibody. Positive IP detected in mouse kidney tissue. Positive WB detected in mouse kidney tissue, HEK-293 cells, human small intestine tissue, L02 cells, mouse small intestine tissue, mouse testis tissue. Positive IHC detected in human colon cancer tissue, human kidney tissue. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of TPK1 recombinant protein (Accession Number: BC014552). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:200-1:2000

IHC: 1:10-1:100

Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116266(TPK1 antibody) at dilution of 1:10



(under 10x lens)

mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:116266(TPK1 Antibody) at dilution of 1:1000

TPK1 antibody

Catalog Number: 116266

IP Result of anti-TPK1 (IP:Catalog No:116266, 4ug; Detection:Catalog No:116266 1:500) with mouse kidney tissue lysate 4000ug.

