

ZSCAN20 antibody



Catalog Number: 117266

Product name

ZSCAN20 antibody

Specificity

Human; other species not tested.

Antibody description

ZSCAN20 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, Jurkat cells. Positive IHC detected in human small intestine tissue.

Observed molecular weight by Western-blot:
100-110 kDa

Preparation

This antibody was obtained by immunization of ZSCAN20 recombinant protein (Accession Number: BC008827). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:117266(ZSCAN20 antibody) at dilution of 1:1500



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:117266(ZSCAN20 Antibody) at dilution of 1:100