

# Anti-SHP1 (2B7) Mouse antibody



Catalog Number: 168266

## Product name

Anti-SHP1 (2B7) Mouse antibody

## Specificity

Human

## Antibody description

SHP1 (2B7) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human SHP-1 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a


## Applications

WB, IP


## Dilutions

WB: 1/1000

## Validations



Western blot detection of SHP-1 in Human tonsil, Daudi, Jurkat and Raji cell lysates using SHP-1 mouse mAb (1:1000 diluted). Predicted band size: 67KDa. Observed band size: 67KDa.



Immunoprecipitation analysis of Raji cell lysates using SHP-1 mouse mAb.