

# ARHGAP27 antibody



Catalog Number: 108266

## Product name

ARHGAP27 antibody

## Specificity

Human; other species not tested.

## Antibody description

ARHGAP27 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in A549 cells. Observed molecular weight by Western-blot: 98 kDa

## Preparation

This antibody was obtained by immunization of ARHGAP27 recombinant protein (Accession Number: XM\_047435551). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:108266(ARHGAP27 Antibody) at dilution of 1:1000



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:108266(ARHGAP27 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)