

# CEP152 antibody



Catalog Number: 109266

## Product name

CEP152 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CEP152 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 110 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_014985). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:2000

## Validations



IP Result of anti-CEP152 (IP:Catalog No:109266, 3ug; Detection:Catalog No:109266 1:500) with



mouse brain tissue lysate 9500ug. mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109266(CEP152 Antibody) at dilution of 1:300