

# CD79B Polyclonal Antibody



Catalog Number: 163666

## Product name

CD79B Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to CD79B

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 29-159 of human CD79B (NP\_000617.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

## Validations



Western blot - CD79B Polyclonal Antibody

Western blot analysis of extracts of Raji cells, using CD79B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .

Immunofluorescence - CD79B Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using CD79B antibody . Blue: DAPI for nuclear staining.