

Anti-FGFR2 antibody

Catalog Number: 176666

Product name

Anti-FGFR2 antibody

Specificity

Human, Mouse

Antibody description

Mouse monoclonal antibody to FGFR2

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, ICC, FC

Dilutions

ICC: 1:50 WB: 1:2000

Validations




Fig1: Western blot analysis on K562 cell lysates using anti-FGFR2 mouse mAb.




Fig2: ICC staining FGFR2 in MCF-7 cells (red). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.