

# A2LD1 antibody

Catalog Number: 107666

---

## Product name

A2LD1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

A2LD1 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, mouse bladder tissue, mouse brain tissue. Observed molecular weight by Western-blot: 17kd

## Preparation

This antibody was obtained by immunization of A2LD1 recombinant protein (Accession Number: XM\_047430708). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:107666(A2LD1 antibody) at dilution of 1:500