

Histone H1.0 antibody

Catalog Number: 111366

Product name

Histone H1.0 antibody

IP: 1:200-1:2000

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Histone H1.0 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IP detected in A431 cells. Positive WB detected in A431 cells, HeLa cells, Jurkat cells, mouse spleen tissue, rat spleen tissue. Observed molecular weight by Western-blot: 32kd

Preparation

This antibody was obtained by immunization of Histone H1.0 recombinant protein (Accession Number: NM_005318). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:111366(H1F0 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells using Catalog No:111366(Histone H1.0 Antibody) at dilution of 1:50 and Rhodamine-Goat anti-Rabbit



IgG

IP Result of anti-Histone H1.0 (IP:Catalog

Histone H1.0 antibody



Catalog Number: 111366

No:111366, 4ug; Detection:Catalog No:111366

1:500) with A431 cells lysate 2400ug.