

Anti-RNF144A antibody



Catalog Number: 176366

Product name

Anti-RNF144A antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to RNF144A

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human rnf144a. aa range:151-200

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1:: Western Blot analysis of various cells using UBCE7IP4 Polyclonal Antibody. Secondary antibody (catalog#: HA1001) was diluted at 1:20000