

TSPAN9 antibody

Catalog Number: 116366

Product name

TSPAN9 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TSPAN9 Rabbit Polyclonal antibody. Positive WB detected in mouse peripheral blood leukocyte tissue, mouse spleen tissue, rat spleen tissue. Positive IP detected in mouse spleen tissue. Observed molecular weight by Western-blot: 29-33 kDa

Preparation

This antibody was obtained by immunization of TSPAN9 recombinant protein (Accession Number: XM_047428126). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:2000

Validations



mouse peripheral blood leukocyte tissue were subjected to SDS PAGE followed by western blot with Catalog No:116366(TSPAN9 Antibody) at dilution of 1:300



IP Result of anti-TSPAN9 (IP:Catalog No:116366, 4ug; Detection:Catalog No:116366 1:500) with mouse spleen tissue lysate 4800ug.