

# IL17RC antibody



Catalog Number: 111766

## Product name

IL17RC antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

IL17RC Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 45-50kd

## Preparation

This antibody was obtained by immunization of IL17RC recombinant protein (Accession Number: BC006411). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111766(IL17RC antibody) at dilution of 1:200