

# Synaptogyrin-4 antibody



Catalog Number: 115766

## Product name

Synaptogyrin-4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Synaptogyrin-4 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HepG2 cells. Observed molecular weight by Western-blot: 31-35kd

## Preparation

This antibody was obtained by immunization of Synaptogyrin-4 recombinant protein (Accession Number: XM\_005258694). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115766(SYNGR4 antibody) at dilution of 1:300