

Anti-LAMP2A Rabbit antibody



Catalog Number: 166766

Product name

Anti-LAMP2A Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

LAMP2A Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human LAMP2A

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC, IP

Dilutions


WB: 1/1000

IHC: 1/20-1/100

FC: 1/200

IP: 1/20

Validations



Western blot detection of LAMP2A in 3T3 cell lysates using LAMP2A Rabbit pAb(1:1000 diluted). Predicted band size:45kDa. Observed band size:110kDa.