

# ADAM28 antibody

Catalog Number: 107766

## Product name

ADAM28 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ADAM28 Rabbit Polyclonal antibody. Positive IP detected in mouse lung tissue. Positive WB detected in mouse lung tissue, BxPC-3 cells, mouse colon tissue, mouse liver tissue, mouse small intestine tissue, mouse spleen tissue, PC-13 cells. Observed molecular weight by Western-blot: 35 kDa

## Preparation

This antibody was obtained by immunization of ADAM28 recombinant protein (Accession Number: NM\_014265). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:107766(ADAM28 antibody) at dilution of 1:800



IP Result of anti-ADAM28 (IP:Catalog No:107766, 4ug; Detection:Catalog No:107766 1:800) with mouse lung tissue lysate 3600ug.