

ASK1 Polyclonal Antibody



Catalog Number: 164076

Product name

ASK1 Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to ASK1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1173-1374 of human ASK1 (NP_005914.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:200

Validations



Immunohistochemistry - ASK1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human stomach using MAP3K5 antibody (40x



lens).

Immunohistochemistry - ASK1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat liver using MAP3K5 antibody (40x lens).



Immunohistochemistry - ASK1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophagus using MAP3K5 antibody (40x lens).