

Anti-Histone H2A.Z Rabbit antibody



Catalog Number: 167076

Product name

Anti-Histone H2A.Z Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

Histone H2A.Z Rabbit polyclonal antibody


Dilutions

WB: 1/1000

Preparation

Antigen: A synthetic peptide of human Histone H2A.Z

Validations



Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Western blot detection of Histone H2A.Z in K562,C6,3T3 cell lysates using Histone H2A.Z Rabbit pAb(1:1000 diluted).Predicted band size:14kDa.Observed band size:14kDa.