

ESR2 Beta-6 antibody



Catalog Number: 110476

Product name

ESR2 Beta-6 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ESR2 Beta-6 Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue. Positive IHC detected in human cervix tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001437). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:110476(ESR2 Beta-6 Antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human cervix using Catalog No:110476(ESR2 Beta-6 antibody) at dilution of 1:50 (under 40x lens)