

MAP3K4 antibody



Catalog Number: 112476

Product name

MAP3K4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MAP3K4 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 170-180 kDa

Preparation

This antibody was obtained by immunization of MAP3K4 recombinant protein (Accession Number: NM_006724). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:112476(MAP3K4 antibody) at dilution of 1:2000



Immunofluorescent analysis of HeLa cells, using MAP3K4 antibody Catalog No:112476 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).