

# Phospho-Tyk2-Y1054/1055 pAb



Catalog Number: 166476

## Product name

Phospho-Tyk2-Y1054/1055 pAb

## Specificity

Human

## Antibody description

Polyclonal antibody to Phospho-Tyk2-Y1054/1055

## Preparation

Antigen: A synthetic phosphorylated peptide around Y1054 & Y1055 of human Tyk2 (NP\_003322.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:1000

## Validations



Western blot - Phospho-Tyk2-Y1054/1055 pAb

Western blot analysis of extracts of HeLa cells, using Phospho-Tyk2-Y1054/1055 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% BSA.