

# Anti-ATG7 Rabbit antibody



Catalog Number: 167476

## Product name

Anti-ATG7 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

ATG7 Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human Atg7(Apg7)

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, ICC/IF

## Dilutions

WB: 1/1000-1/5000

ICC/IF: 1/50-1/200

## Validations



Western blot analysis of Atg7(Apg7) expression in (1) HepG2 cell lysate; (2) Mouse spleen lysate; (3) Rat kidney lysate.