

# Human ZNRF1/NIN283 (His tag) recombinant protein



Catalog Number: 518476

## General Information

### Gene Name Synonym

E3 ubiquitin-protein ligase ZNRF1, nerve injury-induced gene 283 protein, zinc/ring finger protein 1, NIN283, ZNRF1

### Protein Construction

A DNA sequence encoding the Human ZNRF1/NIN283 Pro74-Thr178 is expressed with a 6His tag at the N-terminus, 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL

method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM Tris, 0.15M NaCl, 1mMEDTA, 10% glycerol, pH 8.0 .

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human ZNRF1/NIN283 (His tag) recombinant protein