

# IWS1 antibody

Catalog Number: 111876

## Product name

IWS1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

IWS1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in HeLa cells, human brain tissue, mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 116-120 kDa

## Preparation

This antibody was obtained by immunization of IWS1 recombinant protein (Accession Number: BC017012). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111876(IWS1 Antibody) at dilution of 1:3000



IP Result of anti-IWS1 (IP:Catalog No:111876, 4ug; Detection:Catalog No:111876 1:500) with mouse brain tissue lysate 4000ug.