

Human KLK8/NRPN/PRSS19/TADG14 (His tag) recombinant protein



Catalog Number: 518876

General Information

Gene Name Synonym

Kallikrein-8, hK8, Neuropsin, NP, Ovasin, Serine Protease 19, Serine Protease TADG-14, Tumor-Associated Differentially Expressed Gene 14 Protein, KLK8, NRPN, PRSS19, TADG14

Protein Construction

A DNA sequence encoding the Human KLK8/NRPN/PRSS19/TADG14 Gln29-GLy260 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human KLK8/NRPN/PRSS19/TADG14 (His tag)
recombinant protein