

# JARID2 Polyclonal Antibody



Catalog Number: 166176

## Product name

Rabbit IgG

JARID2 Polyclonal Antibody

## Applications

## Specificity

WB

Human, Mouse, Rat

## Dilutions


WB 1:500 - 1:2000

## Antibody description

## Validations

Polyclonal antibody to JARID2

## Preparation



Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1027-1246 of human JARID2 (NP\_004964.2).

## Formulation

## Western blot - JARID2 Polyclonal Antibody

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Western blot analysis of extracts of various cell lines, using JARID2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type