

# CREST antibody



Catalog Number: 107176

## Product name

CREST antibody

## Immunogen

[Human CREST Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

CREST Mouse Monoclonal antibody. Positive WB detected in MCF-7 cells. Observed molecular weight by Western-blot: 55 kDa

## Preparation

This antibody was obtained by immunization of CREST recombinant protein (Accession Number: NM\_015558). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107176(CREST Antibody) at dilution of 1:500