

# Human B3GAT3 (His tag) recombinant protein



Catalog Number: 519176

## General Information

### Gene Name Synonym

B3GAT3, Galactosylgalactosylxylosylprotein 3-beta-glucuronosyltransferase 3, Beta-1,3-glucuronosyltransferase 3, Glucuronosyltransferase I, GlcAT-I, GlcUAT-I, Gal beta-1,3-Gal-R glucuronosyltransferase,,, REF: C1005

### Protein Construction

A DNA sequence encoding the Human B3GAT3 Glu72-Val335 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, 2mM EDTA, 20% Glycerol, pH 8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human B3GAT3 (His tag) recombinant protein