

# BRCA1 antibody

Catalog Number: 108576

## Product name

BRCA1 antibody

## Immunogen

[Human BRCA1 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

BRCA1 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human lung cancer tissue. Positive WB detected in HeLa cells, A431 cells, HEK-293 cells, MCF-7 cells. Observed molecular weight by Western-blot: 100 kDa

## Preparation

This antibody was obtained by immunization of BRCA1 recombinant protein (Accession Number: BC115037). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117223(BRCA1 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:117223(BRCA1 Antibody) at dilution of 1:50