

Human ACE2/UNQ868/PRO1885 (His tag) recombinant protein



Catalog Number: 518576

General Information

Gene Name Synonym

Angiotensin-Converting Enzyme 2, ACE-Related Carboxypeptidase, Angiotensin-Converting Enzyme Homolog, ACEH, Metalloprotease MPROT15, ACE2

Protein Construction

A DNA sequence encoding the Human ACE2/UNQ868/PRO1885 Gln18-Ser740 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 0.1mM ZnCl₂, 10% Glycerol, pH 7.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human ACE2/UNQ868/PRO1885 (His tag)
recombinant protein