

# Human ACE2/UNQ868/PRO1885 (His tag) recombinant protein



Catalog Number: 518576

## General Information

### Gene Name Synonym

Angiotensin-Converting Enzyme 2, ACE-Related Carboxypeptidase, Angiotensin-Converting Enzyme Homolog, ACEH, Metalloprotease MPROT15, ACE2

### Protein Construction

A DNA sequence encoding the Human ACE2/UNQ868/PRO1885 Gln18-Ser740 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL

method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM TrisHCl, 150mM NaCl, 0.1mM  $\text{ZnCl}_2$ , 10% Glycerol, pH 7.5.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human ACE2/UNQ868/PRO1885 (His tag)  
recombinant protein