

Anti-FBN1 antibody

Catalog Number: 175976

Product name

Anti-FBN1 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to FBN1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human fibrillin-1. aa range:2811-2860

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of mouse-liver cells using FBN1 Polyclonal Antibody