

# Anti-CDK5RAP3 Rabbit antibody



Catalog Number: 166976

## Product name

Anti-CDK5RAP3 Rabbit antibody

WB: 1/1000

IHC: 1/20

## Specificity

Human, Mouse

FC: 1/20

IP: 1/20

## Antibody description

CDK5RAP3 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human CDK5RAP3

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


WB, IHC-P, FC, IP

## Dilutions

## Validations



Immunohistochemistry of CDK5RAP3 in paraffin-embedded Human breast cancer tissue using CDK5RAP3 Rabbit pAb at dilution 1/100



Western blot detection of CDK5RAP3 in K562, Hela cell lysates using CDK5RAP3 Rabbit pAb(1:1000 diluted). Predicted band size: 57kDa. Observed band size: 70kDa.