

HDAC3-specific antibody



Catalog Number: 106976

Product name

HDAC3-specific antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

HDAC3-specific Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, human stomach tissue, U-937 cells. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003883). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:111376(HDAC3-specific antibody) at dilution of 1:200