

ELF2 antibody

Catalog Number: 110276

Product name

ELF2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ELF2 Rabbit Polyclonal antibody. Positive IP detected in A2780 cells. Positive WB detected in A2780 cells, human brain tissue, human liver tissue. Observed molecular weight by Western-blot: 63-64kd, 56-57kd

Preparation

This antibody was obtained by immunization of ELF2 recombinant protein (Accession Number: NM_006874). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:110276(ELF2 antibody) at dilution of 1:500



IP Result of anti-ELF2 (IP:Catalog No:110276, 4ug; Detection:Catalog No:110276 1:500) with A2780 cells lysate 960ug.