

# NRF1/nuclear respiratory factor 1 antibody



Catalog Number: 113276

## Product name

NRF1/nuclear respiratory factor 1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NRF1/nuclear respiratory factor 1 Rabbit Polyclonal antibody. Positive IP detected in mouse skeletal muscle tissue. Positive WB detected in mouse skeletal muscle tissue, human lung tissue. Observed molecular weight by Western-blot: 67 kDa

## Preparation

This antibody was obtained by immunization of NRF1/nuclear respiratory factor 1 recombinant protein (Accession Number: NM\_001293163). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:113276(NRF1 antibody) at dilution of 1:400



IP Result of anti-NRF1 (IP:Catalog No:113276, 4ug; Detection:Catalog No:113276 1:300) with mouse skeletal muscle tissue lysate 3000ug.