

TCF12 Polyclonal Antibody



Catalog Number: 164276

Product name

TCF12 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to TCF12

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 443-682 of human TCF12 (NP_003196.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions


WB 1:200 - 1:1000

Validations



Western blot - TCF12 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TCF12 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



Western blot - TCF12 Polyclonal Antibody

Western blot analysis of extracts of various cells, using TCF12 antibody at 1:1000 dilution ratio through one-step method. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:4000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3min.